

Experiment - 1

Aim: To calculate the number of tyrosine and tryptophan residues in the given BSA samples using UV-Vis spectrophotometer.

Theory: UV-Vis spectroscopy is an analytical technique that measures the amount of discrete wavelengths of UV or visible light are absorbed by or transmitted through a sample in comparison to a reference or blank sample. This property is influenced by the sample composition providing information on the concentration of the sample.

Light has a certain amount of energy which is inversely proportional to its wavelength. A specific amount of energy is needed to promote electrons in a substance to a higher energy state which we can detect as absorption. This is why absorption of light occurs for different wavelengths in different substances.

In this experiment, BSA sample is denatured using 8M urea and its absorbance is measured at 280nm and 283 nm. From the information, the value of extinction coefficient is calculated for both 280nm and 283 nm absorbance values along with the number of tyrosine and tryptophan residues. From the given BSA sequences, the number of tyrosine and tryptophan residue per protein are calculated.

Material Required: BSA sample, 8M urea, 1X PBS (prepared from 10X PBS stock solution by adding 22.5 ml water to 2.5 ml 10X PBS), cuvette, spectrophotometer.

Procedure:

1. At first, 1X PBS is made from 10X PBS stock solution.
2. In spectrophotometer instrument blank is set up with 1X PBS in cuvette to be zero at 280 nm.
3. BSA is now dissolved in PBS.
4. Absorbance is measured at 280 nm with native BSA dissolved in PBS.
5. Now 8M urea (12.016 g of urea in 25 ml 1X PBS) solution is made.
6. Now blank is set to zero with 8M urea solution (urea + 1X PBS).
7. Then BSA is denatured by dissolving in 8M urea solution.
8. Absorbance of the denatured BSA solution is measured at 280 nm.
9. The above procedure is repeated and OD is taken at 283 nm as well.

Observation Table:

Wavelength	PBS (Blank)	PBS + Native BSA	Blank - PBS + Urea	Denatured BSA (+ urea)
A (280 nm)	0	1.011	0	0.699
A (283 nm)	0	0.708	0	0.370

Teacher's Signature _____

Calculations: $c = (1/l) / 66463 = 1.505 \times 10^{-5}$
 $l = 1 \text{ cm} = 10^{-2} \text{ m}$

For 280 nm, $A_1 = \epsilon_1 c l$

$$\epsilon_1 = \frac{A_1}{c l} = \frac{0.699}{1.505 \times 10^{-5} \text{ M} \times 1 \text{ cm}}$$

$$\epsilon_1 = 4.646 \times 10^4 \text{ M}^{-1} \text{ cm}^{-1}$$

For 288 nm,

$$A_2 = \epsilon_2 c l$$

$$\epsilon_2 = \frac{A_2}{c l} = \frac{0.370}{1.505 \times 10^{-5} \text{ M} \times 1 \text{ cm}}$$

$$\epsilon_2 = 2.460 \times 10^4 \text{ M}^{-1} \text{ cm}^{-1}$$

$$\epsilon_1 (280 \text{ nm}) = N_{\text{Trp}} \times 5690 + N_{\text{Tyr}} \times 1280$$

$$4.646 \times 10^4 = N_{\text{Trp}} \times 5690 + N_{\text{Tyr}} \times 1280 \quad \text{--- (1)}$$

$$\epsilon_2 (288 \text{ nm}) = N_{\text{Trp}} \times 4815 + N_{\text{Tyr}} \times 385$$

$$2.46 \times 10^4 = N_{\text{Trp}} \times 4815 + N_{\text{Tyr}} \times 385 \quad \text{--- (2)}$$

$$\textcircled{1} \times 4815 - \textcircled{2} \times 5690$$

$$N_{\text{Tyr}} = \frac{8376.893994 \times 10^4}{3972550} \approx 21$$

Teacher's Signature _____

$$N_{Tyr} \approx 21$$

Putting N_{Tyr} in eqⁿ ①

$$N_{Trp} = \frac{19577.622872}{5600}$$

$$N_{Trp} \approx 3$$

Conclusion:

The values of extinction coefficient are $11.646 \times 10^4 \text{ M}^{-1} \text{ cm}^{-1}$ and $2.46 \times 10^4 \text{ M}^{-1} \text{ cm}^{-1}$ at 280 nm and 298 nm respectively.

The number of tyrosine and tryptophan residues is equal to 21 and 3. Absorbance by protein in denatured form is less than native.

Discussion:

In UV-Vis spectroscopy the wavelength corresponding to the maximum absorbance of the sample is chosen. This ensures maximum sensitivity because the largest response is obtained for a certain analyte concentration.

The technique allows samples to be reused. The instrument is easy to use in which measurements can be made easily.

Teacher's Signature _____