

15/01/2024

Next-Generation Sequencing

- Roche 454 (2006-2013)
- Illumina → most popular (90%)
- PacBio SMRT → single-molecule real-time
- Ion Torrent
- Oxford Nanopore → latest technology
- Data Analysis
- Genomics
- Transcriptomics
- Assembly Algorithms

Sanger sequencing method:

di-deoxy termination method.



capillary electrophoresis

- ✓ one DNA fragment at a time → Drawbacks
- ✓ we get around 1 kb (10^3 bp)
- ✓ took 13 years to sequence the entire human genome!

Human Genome Project

Shotgun sequencing

NGS Technologies:

- ✓ direct sequencing
 - ✓ high-throughput, accurate
 - ✓ cost-effective ~ \$1000 & can go down
- ~ 96 fragments together

Atlas Project

expression of genes changes during diseases

dataset to correlate genetics & drugs

→ The Cancer Genome Atlas (TCGA)
useful in patient-specific therapy.




Precision Medicine

Reads :

- ✓ output sequences from sequencing platform
- ✓ of different lengths

Roche 454 Sequencing :

Pyrosequencing



synthetically -
designed DNA
sequences of
known sequences

primer -
binding
sites

as we do not
know the seq.
of fragment yet

Index sequences/
DNA barcode

to identify each
sample differently

Multiplexing

mixing multiple
samples & sequencing
them together

Pyrosequencing :

relies on synthesizing
the complementary strand



→ A

every time there is
an addition of a
base, a light signal
is emitted

if two same
bases
simultaneously

double the
intensity

fragments




isolation of fragments
↓
clonally amplified.

picotitre plates → 100,000 fragments

put in sequencer

Drawbacks

- ✓ many enzymes & substrates required
- ✓ expensive compared to Illumina, etc.



usually, light intensity increases proportionally with repeated bases.

after 8 repeats

signal is saturated

no more proportional increase

might appear as indels (insertion/deletions)

TTTTTTTTT
homopolymeric stretches
~ 5-20

Illumina - Sequencing by Synthesis (SBS)

Reversible termination

(Sanger has irreversible termination)


modified dNTP

dNTP-ter-Fl

fluorophore

one base at a time

n cycles → n length reads




... process repeated.

29/1/24

Illumina (SBS):

Reversible Termination

dNTP-ter-fl



diff. for the four bases



fl1

fl2

fl3

fl4


detached



flow cell

elongation in opp. dir.

bending to attach to B




Bridge Amplification

after many cycles, we get a cluster of one type of sample

camera

cycles




G

A

T

so on

cycle 2




cluster has same bead length

no. of cycles

limit ~ 250 cycles

nth cycle



in a cluster \rightarrow ~100 molecules that are identical \rightarrow detector reads average signal from a cluster

so after 250 genes/bases, these avg. signals can be confusing & hard to decipher

in the adaptor seq.

we add index sequences

eg: ATTGCC

to identify sequences in a mixed sample

Multiplexing

Single-end (SE) and Paired-end (PE) sequencing:

\rightarrow Illumina Data Processing.

↑ Throughput: high reads / lot of data

150 bp: SE seq.

2x150 bp: PE seq. (like in Illumina video).

Advantages:


- cost-effective (not many enzymes are required).
- polymeric stretch \rightarrow seq. one by one, not all at a time like Roche 454.
- highly accurate.

Drawbacks:


- expensive compared to the newer NGS \rightarrow even a μ l costs in thousands.
- short reads (upto 2x250 bp)

Single Molecule Real Time Sequencing (SMRT): \rightarrow [SBS]

sequencing happens inside ZMWs.



everything else is similar to Illumina \leftarrow



immobilized at the bottom

\rightarrow no amplification step is required.

\rightarrow no cluster formation

\rightarrow no averaging of signals \rightarrow that's why 'single molecule'

\rightarrow more sensitive detection

bottom
as \downarrow of light too large

so incoming DNTPs do not produce signal
signal detected only at illumination

Advantages

- long reads
- avg read length of 2-3 kb, some 10-20 kb
- single-molecule resolution

Drawbacks:

- high error rate ($\sim 14\%$)


Circular Consensus Sequence (CCS) reads

\rightarrow 99% accuracy

1. ATGCCTA TT
2. ATCCCTA TT
3. ATGCCAATT
4. ATGCATATT

take the most freq. at a place

ATGCCTA TT } Consensus sequence



multiple repeats
↓
circular synthesis.

polymerase
at the bottom → always illuminated by light → leads to adverse effect on fidelity

↳ high energy state

short reads but
→ high accuracy

use of Illumina
& SMRT together

→ Hybrid Sequencing

↳ large reads but
low accuracy.